

CD200 Mouse anti-Rat Monoclonal (PE) (MRC OX-2) Antibody - LS-C62679 - LSBio

CatalogID:	LS-C62679
Target:	CD200 molecule
Synonyms:	CD200 Antibody, CD200 antigen Antibody, CD200 molecule Antibody, MOX1 Antibody, MOX2 Antibody, OX-2 membrane glycoprotein Antibody, OX-2 Antibody, MRC Antibody, MRC OX-2 antigen Antibody
Host	CD200 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MRC OX-2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD200 antibody was raised against Rat
Specificity:	Mouse anti-rat OX-2 monoclonal antibody recognizes a monomorphic determinant present on rat thymocytes, brain, follicular dendritic cells in lymphoid organs, vascular endothelium, and at low levels on some smooth muscle and B lymphocytes. The purified bra.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	IHC, Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences