

CD8B / CD8 Beta Mouse anti-Rat Monoclonal (PE) (341) Antibody - LS-C62691 - LSBio

CatalogID:	LS-C62691
Target:	CD8b molecule
Synonyms:	CD8B Antibody, CD8b antigen Antibody, CD8b molecule Antibody, CD8B1 Antibody, p37 Antibody, LY3 Antibody, LYT3 Antibody, Leu2 Antibody
Host	CD8B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	341
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD8B / CD8 Beta antibody was raised against Rat
Antigen Type:	Cells
Immunogen:	CD8B / CD8 Beta antibody was raised against rat/mouse T cell hybrids expressing CD8.
Specificity:	Anti-rat CD8b monoclonal antibody reacts with the beta chain of the CD8 differentiation antigen. CD8b is expressed on most thymocytes and mature T cytotoxic/ suppressor cells (MHC class I restricted). While the CD8a and CD8b form a heterodimer on the surf.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Western blot, Immunoprecipitation, Flow Cytometry (1 µg/10E6 cells), Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences