

THY1 / CD90 Mouse anti-Rat Monoclonal (PE) (MRC OX-7) Antibody - LS-C62714 - LSBio	
CatalogID:	LS-C62714
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MRC OX-7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	THY1 / CD90 antibody was raised against Rat
Antigen Type:	Cells
Immunogen:	THY1 / CD90 antibody was raised against rat thymocyte Thy-1 antigen.
Specificity:	Anti-rat CD90 (Thy 1) monoclonal antibody recognizes the Thy 1.1 antigenic determinant, designated CD90 on rat as well as mouse cells. This particular determinant has been defined to be monomorphic within rats but polymorphic in the mouse. The antibody re.
Reactivity:	Rat, Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	