

CD4 Rat anti-Mouse Monoclonal (PE) (YTS 191.1.2) Antibody - LS-C62732 - LSBio

CatalogID:	LS-C62732
Target:	CD4 molecule
Synonyms:	CD4 Antibody, CD4 antigen Antibody, CD4 receptor Antibody, CD4mut Antibody, CD4 antigen (p55) Antibody, CD4 molecule Antibody
Host	CD4 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	YTS 191.1.2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD4 antibody was raised against Mouse
Immunogen:	CD4 antibody was raised against murine thymocytes.
Specificity:	The CD4 (L3/T4) antigen is thought to be the murine equivalent of the human Leu3/OKT4 antigen. This T cell surface molecule appears to be expressed by the helper/inducer subset of murine T cells and by delayed hypersensitivity T cells but not by cytotoxic.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences