

CD45RC Mouse anti-Rat Monoclonal (FITC) (MRC OX-22) Antibody - LS-C62742 - LSBio

CatalogID:	LS-C62742
Target:	CD45RC
Host	CD45RC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MRC OX-22
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45RC antibody was raised against Rat
Immunogen:	CD45RC antibody was raised against phytohemagglutinin-activated rat lymph node cells.
Specificity:	Anti-rat CD45RC monoclonal antibody reacts with the high molecular weight form (190-240 kD) of the leukocyte common antigen designated as CD45RC. This designation includes the 240 kD form on B cells as well as a subfraction of the 190, 200 and 220 kD f.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg or 500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences