

C3AR / C3a Receptor Mouse anti-Rat Monoclonal (PE) (74/rC3aRZ1) Antibody - LS-C62771 - LSBio

CatalogID:	LS-C62771
Target:	complement component 3a receptor 1 (C3AR1)
Synonyms:	C3AR1 Antibody, AZ3B Antibody, C3AR Antibody, C3R1 Antibody, C3a receptor Antibody, C3a-R Antibody, HNFAG09 Antibody, Anaphylatoxin c3a receptor Antibody, C3a anaphylatoxin receptor Antibody
Family / Subfamily:	GPCR / Anaphylatoxin
Host	C3AR1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	74/rC3aRZ1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	C3AR / C3a Receptor antibody was raised against Rat
Specificity:	C3aR binds to anaphylatoxin C3a which, among other things, is involved in mast cell degranulation and recruitment of immune cells to the site of inflammation. Activation of the complement cascade produces various fragments, including C3a and C5a. C3aR is
Reactivity:	Rat
Uses:	Undetermined
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences