

THY1 / CD90 Rabbit anti-Mouse Polyclonal Antibody - LS-C62788 - LSBio

CatalogID:	LS-C62788
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Reactivity:	Mouse
Uses:	Undetermined
Size:	5 ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences