

Mouse IgG + IgA + IgM Rabbit anti-Mouse Polyclonal (FITC) Antibody - LS-C62812 - LSBio

CatalogID:	LS-C62812
Target:	Mouse IgA/IgG/IgM
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse
Specificity:	The antiserum is prepared by immunizing rabbits with mouse IgG. It is polyvalent with respect to classes of mouse immunoglobulin. The antiserum is specific for mouse immunoglobulin: It does not react with mouse serum proteins other than immunoglobulins; i
Reactivity:	Mouse
Purification:	Unpurified serum
Reconstitution:	Distilled Water. Possible additional volumes for resuspension: 1 ml
Presentation:	Lyophilized.
Uses:	Undetermined
Size:	1 ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences