

Rabbit IgG + IgA + IgM Porcine anti-Rabbit Polyclonal (FITC) Antibody - LS-C62816 - LSBio	
<b>CatalogID:</b>	LS-C62816
<b>Target:</b>	Rabbit IgA + IgG + IgM
<b>Host</b>	Porcine
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Rabbit
<b>Specificity:</b>	The antiserum is prepared by immunizing swine with rabbit IgG. The antiserum is specific for rabbit immunoglobulin: It does not react with rabbit serum proteins other than immunoglobulins: it does not cross-react with immunoglobulins of other species
<b>Reactivity:</b>	Rabbit
<b>Purification:</b>	Unpurified serum
<b>Reconstitution:</b>	Distilled Water. Possible additional volumes for resuspension: 1 ml
<b>Presentation:</b>	Lyophilized, Carbonate buffer.
<b>Recommended Storage:</b>	Store at -20°C or below before reconstitution. Aliquot and store at -70° after reconstitution. Prepare aliquots in appropriate volume for single use to avoid repeated freeze-thawing.
<b>Uses:</b>	Undetermined
<b>Size:</b>	1 ml
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 10/2/2014 © 2014 LifeSpan BioSciences	