

TCR D011.10 Mouse anti-Mouse Monoclonal (PE, Cy5) (KJ1-26) Antibody - LS-C62830 - LSBio	
CatalogID:	LS-C62830
Target:	T Cell Receptor D011.10
Host	T Cell Receptor D011.10 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	KJ1-26
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	TCR D011.10 antibody was raised against Mouse
Immunogen:	TCR D011.10 antibody was raised against t cell hybridoma, D011.10.
Specificity:	This monoclonal antibody binds the T-cell receptor (TCR) expressed on DO11.10 and DO1 1.10.24 murine T cell 1 hybridomas, but not on any other T cell hybridomas tested. The TCR expressed on DO11.10 recognizes an ovalbumin peptide in the context of I-Ad ma
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	