

CD55 Hamster anti-Mouse Monoclonal (PE) (RIKO-3) Antibody - LS-C62833 - LSBio

CatalogID:	LS-C62833
Target:	CD55 molecule, decay accelerating factor for complement (Cromer blood group)
Synonyms:	CD55 Antibody, CD55 antigen Antibody, Decay-accelerating factor Antibody, TC Antibody, CR Antibody, CROM Antibody, DAF Antibody
Host	CD55 antibody was produced in Hamster
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	RIKO-3
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD55 antibody was raised against Mouse
Immunogen:	CD55 antibody was raised against gPI-DAF transfected CHO cells.
Specificity:	Anti-Mouse DAF monoclonal antibody (Clone: RIKO-3) detects DAF; one of several membrane proteins that prevent host cells from homologous complement attack. The two mouse DAF isoforms are effective against rat complement, but not against human and guinea p.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Store at +4°C. Do not freeze. Avoid prolonged exposure to light
Uses:	Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences