

TLR2 Mouse anti-Human Monoclonal (PE) (11G7) Antibody - LS-C62846 - LSBio

CatalogID:	LS-C62846
Target:	toll-like receptor 2 (TLR2)
Synonyms:	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	11G7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TLR2 antibody was raised against Human
Immunogen:	TLR2 antibody was raised against human TLR2.
Specificity:	Anti-human Toll-like Receptor 2 (TLR2) monoclonal antibody detects human TLR2 (CD282). This single membrane spanning, 784 amino acid, 90 KDa protein is a type of pattern recognition receptor (PRR) and recognizes structurally conserved microbial molecules
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Immunoprecipitation, Flow Cytometry (1 µg/10E6 cells), ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences