

MHC Class II I-Ab and I-Ad Mouse anti-Mouse Monoclonal (PE) (25-9-17s) Antibody - LS-C62854 - LSBio

CatalogID:	LS-C62854
Target:	MHC Class II I-Ab and I-Ad
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	25-9-17s
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Mouse
Antigen Type:	Cells
Immunogen:	C3H. SW splenocytes.
Specificity:	is a monoclonal antibody that reacts with I-Ab and I-Ad antigens. Cross reaction with H-2p and H-2q was also found
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Store at 4°C. Stable for one year.
Uses:	Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences