

CD45RB Rat anti-Mouse Monoclonal (FITC) (16A) Antibody - LS-C62865 - LSBio	
CatalogID:	LS-C62865
Target:	CD45RB
Host	CD45RB antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	16A
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45RB antibody was raised against Mouse
Immunogen:	CD45RB antibody was raised against th2 cell clones, final boost with TH2 clone D10.
Specificity:	Anti-mouse CD45RB monoclonal antibody reacts with the CD45 isoform containing the exon B dependent epitope. CD45RB is highly expressed on peripheral B cells, cytotoxic T cells, a subset of T helper cells and most thymocytes
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.1% sodium azide, EIA grade BSA.
Uses:	IHC - Paraffin, IHC - Frozen, Immunoprecipitation, Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg or 300 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	