

CD45RB Rat anti-Mouse Monoclonal (PE) (16A) Antibody - LS-C62866 - LSBio

CatalogID:	LS-C62866
Target:	CD45RB
Host	CD45RB antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	16A
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45RB antibody was raised against Mouse
Specificity:	Anti-mouse CD45RB monoclonal antibody reacts with the CD45 isoform containing the exon B dependent epitope. CD45RB is highly expressed on peripheral B cells, cytotoxic T cells, a subset of T helper cells and most thymocytes
Reactivity:	Mouse
Uses:	Undetermined
Size:	50 µg or 300 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences