

Ly6A/E / Sca-1 Rat anti-Mouse Monoclonal (PE) (D7) Antibody - LS-C62868 - LSBio

CatalogID:	LS-C62868
Target:	Ly 6A/E (Sca-1)
Host	Ly 6A/E (Sca-1) antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	D7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Ly6A/E / Sca-1 antibody was raised against Mouse
Specificity:	This monoclonal antibody reacts with Ly-6A.2 and Ly-6E.1, members of the Ly-6 multigene family. Ly-6A/E molecules are expressed on multipotent hematopoietic stem cells1 activated T cells, and peripheral B cells on mice that express the Ly-6.2 haplotype
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, EIA grade BSA.
Uses:	Western blot, Immunoprecipitation, Flow Cytometry (0.25 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 300 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences