

CD59 Mouse anti-Human Monoclonal (PE) (1F5) Antibody - LS-C62882 - LSBio

CatalogID:	LS-C62882
Target:	CD59 molecule, complement regulatory protein
Synonyms:	CD59 Antibody, 16.3A5 Antibody, 1F5 Antibody, 1F5 antigen Antibody, CD59 glycoprotein Antibody, EL32 Antibody, EJ16 Antibody, G344 Antibody, HRF-20 Antibody, MAC-IP Antibody, MACIF Antibody, Ly-6-like protein Antibody, MEM43 Antibody, MIRL Antibody, MEM43 antigen Antibody, MIN2 Antibody, MIN3 Antibody, MSK21 Antibody, p18-20 Antibody, Protectin Antibody, MAC-inhibitory protein Antibody, MIC11 Antibody, T cell-activating protein Antibody, CD59 antigen Antibody, EJ30 Antibody, HRF20 Antibody, Human leukocyte antigen MIC11 Antibody, Lymphocytic antigen CD59/MEM43 Antibody, MIN1 Antibody
Host	CD59 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	1F5
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD59 antibody was raised against Human
Specificity:	CD59 (also known as HRF20, protectin) is a 20 kD glycoprotein attached to the plasma membrane by a glycosylphosphatidylinositol (GPI) anchor. It blocks MAC formation in the complement pathway by interfering with C9 membrane insertion and polymerization.
Reactivity:	Human
Uses:	Undetermined
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences