

**THY1 / CD90 Mouse anti-Mouse Monoclonal (5a-8) Antibody - LS-C62888 - LSBio**

<b>CatalogID:</b>	LS-C62888
<b>Target:</b>	Thy-1 cell surface antigen (THY1)
<b>Synonyms:</b>	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
<b>Host</b>	THY1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b
<b>Clone Name:</b>	5a-8
<b>Immunogen Species:</b>	THY1 / CD90 antibody was raised against Mouse
<b>Immunogen:</b>	THY1 / CD90 antibody was raised against cBA/J.
<b>Specificity:</b>	Anti-mouse CD90 (Thy 1.2) monoclonal antibody reacts with all T lymphocytes from mouse strains expressing the Thy 1.2 phenotype (e. g. C57BL/6, C3H/He, DBA/2, CBA/J, BALB/c), but does react with lymphocytes expressing the Thy 1.1 phenotype (e. g. AKR/J, B6.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Tissue culture supernatant
<b>Reconstitution:</b>	Distilled Water. Possible additional volumes for resuspension: 1 ml
<b>Presentation:</b>	Lyophilized.
<b>Uses:</b>	Cytotoxic Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 ml
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences