

THY1 / CD90 Mouse anti-Mouse Monoclonal (APC) (5a-8) Antibody - LS-C62889 - LSBio	
CatalogID:	LS-C62889
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	5a-8
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Immunogen:	THY1 / CD90 antibody was raised against cBA/J.
Specificity:	Anti-mouse CD90 (Thy 1.2) is a monoclonal antibody that reacts with all T lymphocytes from mouse strains expressing the Thy 1.2 phenotype (e. g. C57BL/6, C3H/He, DBA/2, CBA/J, BALB/c, BXSB, NOD, NZB and SJL), but does not react with lymphocytes expressing.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (1 µg/10E6 cells), Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014 © 2014 LifeSpan BioSciences	