

**MHC Class II I-Ab (public; Ia. m8) Mouse anti-Mouse Monoclonal (PE) (28-16-8s) Antibody - LS-C62896 - LSBio**

<b>CatalogID:</b>	LS-C62896
<b>Target:</b>	MHC Class II I-Ab (public; Ia. m8)
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgM
<b>Clone Name:</b>	28-16-8s
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	Mouse
<b>Specificity:</b>	Anti- I-Ab monoclonal antibody reacts with the I-Ab encoded MHC class II antigen expressed on mouse strains of the H-2b haplotype. It also reacts with the I-Ad encoded MHC class II antigen expressed on mouse strains of the H- 2d haplotype. The reaction pa
<b>Reactivity:</b>	Mouse
<b>Uses:</b>	Undetermined
<b>Size:</b>	50 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences