

RPE Rabbit anti-Gracilaria Polyclonal (FITC) Antibody - LS-C62913 - LSBio	
<b>CatalogID:</b>	LS-C62913
<b>Target:</b>	R-Phycoerythrin
<b>Host</b>	R-Phycoerythrin antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	RPE antibody was raised against Gracilaria
<b>Reactivity:</b>	Gracilaria
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	