

CD5 Rat anti-Mouse Monoclonal (FITC) (53-7.8) Antibody - LS-C62915 - LSBio	
<b>CatalogID:</b>	LS-C62915
<b>Target:</b>	CD5 molecule
<b>Synonyms:</b>	CD5 Antibody, CD5 antigen Antibody, CD5 antigen (p56-62) Antibody, CD5 molecule Antibody, Lymphocyte antigen T1/Leu-1 Antibody, LEU1 Antibody, T1 Antibody
<b>Host</b>	CD5 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	53-7.8
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD5 antibody was raised against Mouse
<b>Reactivity:</b>	Mouse
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	