

**CD5 Rat anti-Mouse Monoclonal (PE) (53-7.3) Antibody - LS-C62916 - LSBio**

<b>CatalogID:</b>	LS-C62916
<b>Target:</b>	CD5 molecule
<b>Synonyms:</b>	CD5 Antibody, CD5 antigen Antibody, CD5 antigen (p56-62) Antibody, CD5 molecule Antibody, Lymphocyte antigen T1/Leu-1 Antibody, LEU1 Antibody, T1 Antibody
<b>Host</b>	CD5 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	53-7.3
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD5 antibody was raised against Mouse
<b>Reactivity:</b>	Mouse
<b>Uses:</b>	Undetermined
<b>Size:</b>	50 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences