

TNFSF4 / OX40L Rat anti-Mouse Monoclonal (FITC) (OX-89) Antibody - LS-C62948 - LSBio	
CatalogID:	LS-C62948
Target:	tumor necrosis factor (ligand) superfamily, member 4 (TNFSF4)
Synonyms:	TNFSF4 Antibody, CD134 ligand Antibody, CD134L Antibody, CD252 antigen Antibody, Glycoprotein Gp34 Antibody, OX-40L Antibody, OX40 ligand Antibody, TXGP1 Antibody, OX40L Antibody, OX4OL Antibody, CD252 Antibody, gp34 Antibody, OX40 antigen ligand Antibody
Family / Subfamily:	TNF / not assigned-TNF
Host	TNFSF4 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	OX-89
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TNFSF4 / OX40L antibody was raised against Mouse
Immunogen:	TNFSF4 / OX40L antibody was raised against cD4. CD134L chimeric protein.
Specificity:	Anti-mouse CD134L (OX40L) monoclonal antibody reacts with a type II membrane protein that is similar to TNF that is found on activated antigen presenting cells (APC's) including B cells, dendritic cells and macrophages. CD134L mRNA is expressed in various.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze-thaw cycles and prolonged exposure to light.
Uses:	IHC - Frozen, Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	