

Basigin / Emmprin / CD147 Rat anti-Mouse Monoclonal (PE) (OX114) Antibody - LS-C62953 - LSBio

CatalogID:	LS-C62953
Target:	basigin (Ok blood group) (BSG)
Synonyms:	BSG Antibody, 5F7 Antibody, Basigin Antibody, Basigin (Ok blood group) Antibody, CD147 antigen Antibody, Collagenase stimulatory factor Antibody, EMMPRIN Antibody, OK Antibody, OK blood group antigen Antibody, TCSF Antibody, CD147 Antibody, M6 Antibody
Host	BSG antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	OX114
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Basigin / Emmprin / CD147 antibody was raised against Mouse
Immunogen:	Basigin / Emmprin / CD147 antibody was raised against recombinant md0-2rCD4d3+4.
Specificity:	CD147 is a widely expressed membrane glycoprotein (also known as OX47, basigin, EMMPRIN and HT7; expressed on leukocytes, red cells, platelets and endothelial cells) and has been implicated in matrix metalloproteinase induction, cell adhesion, retinal cell.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences