

PDCD1 / CD279 / PD-1 Hamster anti-Mouse Polyclonal (PE) Antibody - LS-C62961 - LSBio	
<b>CatalogID:</b>	LS-C62961
<b>Target:</b>	programmed cell death 1 (PDCD1)
<b>Synonyms:</b>	PDCD1 Antibody, CD279 Antibody, HPD-1 Antibody, HPD-I Antibody, PD1 Antibody, Protein PD-1 Antibody, SLEB2 Antibody, PD-1 Antibody, CD279 antigen Antibody, Programmed cell death 1 Antibody
<b>Family / Subfamily:</b>	Apoptosis
<b>Host</b>	PDCD1 antibody was produced in Hamster
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	PDCD1 / CD279 / PD-1 antibody was raised against Mouse
<b>Immunogen:</b>	PDCD1 / CD279 / PD-1 antibody was raised against pD-1 cDNA transfectant.
<b>Specificity:</b>	PD1 monoclonal antibody reacts with mouse PD-1 (programmed death-1), a 55 kD member of the Ig superfamily that has been demonstrated to play a key role in the regulation of immune response and peripheral tolerance. PD1 is expressed on activated T and B I.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.02% sodium azide, EIA grade BSA.
<b>Uses:</b>	Flow Cytometry (2 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	