

Mouse IgG Goat anti-Mouse Polyclonal (FITC) Antibody - LS-C63003 - LSBio	
<b>CatalogID:</b>	LS-C63003
<b>Target:</b>	Mouse IgG
<b>Host</b>	Mouse IgG antibody was produced in Goat
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Mouse IgG antibody was raised against Mouse
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	10 mM sodium phosphate, 0.15 M sodium chloride, 0.05% sodium azide, EIA grade BSA
<b>Recommended Storage:</b>	+4°C or -20°C, Avoid repeated freezing and thawing.
<b>Uses:</b>	Flow Cytometry (0.25 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	2 ml
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	