

Mouse IgG Goat anti-Mouse Polyclonal (FITC) Antibody - LS-C63008 - LSBio	
CatalogID:	LS-C63008
Target:	Mouse IgG
Host	Mouse IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse IgG antibody was raised against Mouse
Reactivity:	Mouse
Purification:	Immunoaffinity purified
Presentation:	0.1% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	1 ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	