

Mouse IgG Goat anti-Mouse Polyclonal (PE) Antibody - LS-C63009 - LSBio

CatalogID:	LS-C63009
Target:	Mouse IgG
Host	Mouse IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Mouse IgG antibody was raised against Mouse
Reactivity:	Mouse
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide, EIA grade BSA.
Recommended Storage:	Store at +4°C. DO NOT FREEZE.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	0.5 ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences