

Rat IgG Goat anti-Rat Polyclonal Antibody - LS-C63012 - LSBio

CatalogID:	LS-C63012
Target:	Rat IgG
Host	Rat IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	Rat IgG antibody was raised against Rat
Reactivity:	Rat
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences