

Rat IgG Goat anti-Rat Polyclonal (APC) Antibody - LS-C63015 - LSBio	
CatalogID:	LS-C63015
Target:	Rat IgG
Host	Rat IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	Rat IgG antibody was raised against Rat
Reactivity:	Rat
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide, EIA grade BSA.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	0.5 ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	