

Rat IgG Goat anti-Rat Polyclonal (HRP) Antibody - LS-C63016 - LSBio	
CatalogID:	LS-C63016
Target:	Rat IgG
Host	Rat IgG antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Immunogen Species:	Rat IgG antibody was raised against Rat
Reactivity:	Rat
Presentation:	PBS, 40% glycerol, EIA grade BSA.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	1 ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	