

**Rat IgG Goat anti-Rat Polyclonal (Biotin) Antibody - LS-C63017 - LSBio**

<b>CatalogID:</b>	LS-C63017
<b>Target:</b>	Rat IgG
<b>Host</b>	Rat IgG antibody was produced in Goat
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Biotin
<b>Immunogen Species:</b>	Rat IgG antibody was raised against Rat
<b>Reactivity:</b>	Rat
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, EIA grade BSA.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1 ml
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences