

**HLA-A/B/C Mouse anti-Human Monoclonal (PE) (W6/32 HLK) Antibody - LS-C63035 - LSBio**

<b>CatalogID:</b>	LS-C63035
<b>Target:</b>	MHC Class I HLA-A/HLA-B/HLA-C
<b>Host</b>	MHC Class I HLA-A/HLA-B/HLA-C antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	W6/32 HLK
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	HLA-A/B/C antibody was raised against Human
<b>Specificity:</b>	Anti-human HLA-ABC monoclonal antibody recognizes an epitope common among 43 kDa chains of the HLA-ABC antigens. These antigens appear on virtually every human nucleated cell. T
<b>Reactivity:</b>	Human
<b>Uses:</b>	Undetermined
<b>Size:</b>	50 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences