

**HLA-DR Rat anti-Human Monoclonal (FITC) (YD1/63.4.10) Antibody - LS-C63041 - LSBio**

<b>CatalogID:</b>	LS-C63041
<b>Target:</b>	MHC Class II HLA-DR
<b>Host</b>	MHC Class II HLA-DR antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	YD1/63.4.10
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	HLA-DR antibody was raised against Human
<b>Immunogen:</b>	HLA-DR antibody was raised against dAUDI cells.
<b>Specificity:</b>	FITC anti-human HLA-DR monoclonal antibody recognizes the HLADR (MHC class II) antigen
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.02% sodium azide, EIA grade BSA.
<b>Uses:</b>	Flow Cytometry (0.1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences