

HLA-DR Rat anti-Human Monoclonal (PE) (YD1/63.4.10) Antibody - LS-C63042 - LSBio

CatalogID:	LS-C63042
Target:	MHC Class II HLA-DR
Host	MHC Class II HLA-DR antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	YD1/63.4.10
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HLA-DR antibody was raised against Human
Specificity:	FITC anti-human HLA-DR monoclonal antibody recognizes the HLADR (MHC class II) antigen
Reactivity:	Human
Uses:	Undetermined
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences