

HMOX1 / HO-1 Mouse anti-Rat Monoclonal (PE) (HO-1-2) Antibody - LS-C63206 - LSBio	
CatalogID:	LS-C63206
Target:	heme oxygenase (decycling) 1 (HMOX1)
Synonyms:	HMOX1 Antibody, BK286B10 Antibody, HO-1 Antibody, HSP32 Antibody, Heat shock protein, 32-kD Antibody, Heme oxygenase (decycling) 1 Antibody, Heme oxygenase 1 Antibody, HO1 Antibody
Host	HMOX1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	HO-1-2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HMOX1 / HO-1 antibody was raised against Rat
Immunogen:	HMOX1 / HO-1 antibody was raised against native rat HO-1 protein.
Reactivity:	Rat, Human, Monkey, Mouse, Dog, Guinea pig, Hamster, Pig
Purification:	Protein G purified
Presentation:	PBS, pH 7.5, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (1:100) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	