

**HSPA8 / HSC70 Rat anti-Hamster Monoclonal (PE) (1B5) Antibody - LS-C63254 - LSBio**

<b>CatalogID:</b>	LS-C63254
<b>Target:</b>	heat shock 70kDa protein 8 (HSPA8)
<b>Synonyms:</b>	HSPA8 Antibody, Heat shock 70 kDa protein 8 Antibody, Heat shock 70kd protein 10 Antibody, Heat shock 70kD protein 8 Antibody, Heat shock cognate protein 54 Antibody, HSP71 Antibody, HSPA10 Antibody, HSC71 Antibody, HSP73 Antibody, LPS-associated protein 1 Antibody, Heat shock 70kDa protein 8 Antibody, NIP71 Antibody, HSC54 Antibody, HSC70 Antibody, LAP1 Antibody
<b>Host</b>	HSPA8 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	1B5
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	HSPA8 / HSC70 antibody was raised against Hamster
<b>Immunogen:</b>	HSPA8 / HSC70 antibody was raised against native hamster Hsc70 (Hsp73) protein.
<b>Reactivity:</b>	Hamster, Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Pig, Rabbit, Sheep
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% azide
<b>Recommended Storage:</b>	Store at 4°C.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences