

**HSP90 / Heat Shock Protein 90 Mouse anti-Achlya ambisexualis Monoclonal (PE) (AC88) Antibody
 - LS-C63265 - LSBio**

CatalogID:	LS-C63265
Target:	HSP90 / Heat Shock Protein 90
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	AC88
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Achlya ambisexualis
Immunogen:	Achlya ambisexualis (water mold) Hsp90 protein.
Reactivity:	Achlya ambisexualis, Human, Monkey, Mouse, Rat, Dog, Guinea pig, Hamster, Pig, Rabbit, Sheep, Fish, Bacteria, Beluga, C. elegans, Mold, Scallop
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (1:100) (Optimal dilution to be determined by the researcher)
Size:	50 µg or 200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences