

FITC Rabbit Polyclonal Antibody - LS-C66242 - LSBio

CatalogID:	LS-C66242
Target:	Fluorescein Isothiocyanate (FITC)
Host:	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Antigen Species:	
Antigen Type:	Hapten
Immunogen:	Fluorescein isothiocyanate-KLH.
Specificity:	Recognizes fluorescein isothiocyanate: 100%.
Purification:	Antiserum
Presentation:	Sterile-filtered, defibrinated neat serum. No preservatives added.
Usage Summary:	Suitable for use in ELISA and Immunoprecipitation. Immunoprecipitation: Binds ~8 ug FITC/ml.
Uses:	Immunoprecipitation (IP), ELISA (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 3/15/2012

© 2012 LifeSpan BioSciences