

Soybean Agglutinin Rabbit anti-Soybean Polyclonal (FITC) Antibody - LS-C66788 - LSBio	
CatalogID:	LS-C66788
Target:	Soybean Agglutinin (SBA)
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Soybean
Antigen Type:	Purified protein
Immunogen:	Purified Glycine max agglutinin (SBA, soybean agglutinin).
Specificity:	Recognizes soybean agglutinin.
Reactivity:	Soybean
Purification:	Ammonium sulfate precipitation
Presentation:	PBS, pH 7.3, sodium azide, 40% glycerol.
Usage Summary:	Suitable for use in Immunodiffusion and IEP.
Uses:	Immunodiffusion, Immunoelectrophoresis (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	