

Canine Heartworm (Dirofilaria immitis) Goat anti-Dirofilaria immitis Polyclonal (FITC) Antibody - LS-C67259 - LSBio

CatalogID:	LS-C67259
Target:	Canine Heartworm (Dirofilaria immitis)
Host	Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Dirofilaria immitis
Antigen Type:	Cells
Immunogen:	Dirofilaria immitis.
Specificity:	Soluble alpha structural antigens of the adult worm. Antiserum is not absorbed for related filarial Nematodes. Uninfected Cell Reactivity: Non-reactive with canine serum.
Reactivity:	Dirofilaria immitis
Purification:	Purified
Presentation:	PBS, pH 7.2, BSA (10 mg/ml), 0.1% sodium azide.
Usage Summary:	Suitable for use in ELISA and Immunofluorescence. Immunofluorescence: 1:10-1:50. Direct FA staining of target antigens in a permissive tissue culture system. Acetone fixation of the antigen source is recommended prior to staining.
Uses:	Immunofluorescence (1:10 - 1:50), ELISA (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences