

Ia Antigen Mouse anti-Chicken Monoclonal (PE) Antibody - LS-C68318 - LSBio

CatalogID:	LS-C68318
Target:	Ia Antigen
Host	Ia Antigen antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgM,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Ia Antigen antibody was raised against Chicken
Specificity:	Recognizes a 65kD glycoprotein consisting of a 32kD alpha-chain non-covalently bound to a 27kD beta-chain.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Identification and enumeration of Ia ⁺ cells by flow cytometry. Immunoprecipitation. Inhibition of MLR. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences