

ALB / Serum Albumin Goat anti-Human Polyclonal (FITC) Antibody - LS-C68850 - LSBio

CatalogID:	LS-C68850
Target:	albumin
Synonyms:	ALB Antibody, Albumin (32 AA) Antibody, Albumin (AA 34) Antibody, Albumin Antibody, Growth-inhibiting protein 20 Antibody, PRO1341 Antibody, PRO0883 Antibody, Serum albumin Antibody, PRO0903 Antibody
Host	ALB antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ALB / Serum Albumin antibody was raised against Human
Specificity:	Recognizes human Albumin Serum.
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.4, 0.2% BSA, 0.05% sodium azide.
Recommended Storage:	+4°C (short-term), -20°C (long-term). Protect from light.
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.5 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences