

PNA / Peanut Agglutinin Rabbit anti-Peanut Polyclonal (FITC) Antibody - LS-C69216 - LSBio

CatalogID:	LS-C69216
Target:	PNA / Peanut Agglutinin
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Peanut
Antigen Type:	Protein
Immunogen:	Arachis hypogaea (PNA, peanut).
Specificity:	Recognizes Arachis hypogaea (PNA, peanut).
Reactivity:	Peanut
Purification:	Ammonium sulfate precipitation
Presentation:	PBS, pH 7.2, 0.05% sodium azide
Recommended Storage:	Short term 4 C, Long term -20 C, Product is photosensitive and should be protected from light. Avoid freeze-thaw cycles.
Usage Summary:	Suitable for use in Immunodiffusion, Immunohistochemistry and IEP.
Uses:	IHC, Immunodiffusion, Immunoelectrophoresis (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Concentration:	1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences