

Mouse IgG2a Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C70282 - LSBio

CatalogID:	LS-C70282
Target:	Mouse IgG2a
Host	Mouse IgG2a antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse IgG2a antibody was raised against Mouse
Antigen Type:	Purified protein
Immunogen:	Mouse IgG2a antibody was raised against purified mouse IgG2a.
Specificity:	Recognizes mouse gamma 2a heavy chain, binding to all allotypes.
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.1% sodium azide, 50% glycerol.
Usage Summary:	Suitable for use in Immunomicroscopy and Flow Cytometry and FACS Analysis. Flow Cytometry: 5 ug/ml.
Uses:	Immunofluorescence, Flow Cytometry (5 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences