

**Rat IgG2a Rat anti-Rat Monoclonal (HRP) Antibody - LS-C70283 - LSBio**

<b>CatalogID:</b>	LS-C70283
<b>Target:</b>	Rat IgG2a
<b>Host</b>	Rat IgG2a antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Conjugations:</b>	Horseradish Peroxidase (HRP)
<b>Immunogen Species:</b>	Rat IgG2a antibody was raised against Rat
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Rat IgG2a antibody was raised against purified Rat IgG2a heavy chain.
<b>Specificity:</b>	Recognizes mouse gamma 2a immunoglobulin heavy chain and does not cross-react with other mouse immunoglobulin classes or subclasses.
<b>Reactivity:</b>	Rat, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, 0.01% Thimerosal, 50% glycerol.
<b>Usage Summary:</b>	ELISA: 500 ng/ml.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences