

IL4R / CD124 Goat anti-Mouse Polyclonal (Extracellular Domain) (PE) Antibody - LS-C70897 - LSBio

CatalogID:	LS-C70897
Target:	interleukin 4 receptor (IL4R)
Synonyms:	IL4R Antibody, CD124 Antibody, CD124 antigen Antibody, IL-4RA Antibody, IL4R nirs variant 1 Antibody, IL4RA Antibody, Interleukin 4 receptor Antibody, IL-4 receptor subunit alpha Antibody, IL-4R subunit alpha Antibody, IL-4R-alpha Antibody
Family / Subfamily:	Interleukin
Host	IL4R antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	IL4R / CD124 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	IL4R / CD124 antibody was raised against nS0-derived rmlL-4R extracellular domain.
Specificity:	Recognizes mouse IL-4R.
Epitope:	Extracellular Domain
Reactivity:	Mouse
Purification:	Immunoaffinity purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: Neat. Use 10 ul to label 1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences