

IL22 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C71010 - LSBio

CatalogID:	LS-C71010
Target:	interleukin 22 (IL22)
Synonyms:	IL22 Antibody, Cytokine Zcyto18 Antibody, ILTIF Antibody, Interleukin 22 Antibody, IL-22 Antibody, IL-D110 Antibody, TIFIL-23 Antibody, Zcyto18 Antibody, IL-TIF Antibody, Interleukin-22 Antibody
Family / Subfamily:	Interleukin
Host	IL22 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	IL22 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	IL22 antibody was raised against e. coli-derived rmlL-22.
Specificity:	Recognizes mouse IL-22.
Reactivity:	Mouse
Predicted Reactivity:	Human, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.5% BSA, 0.1% sodium azide.
Recommended Storage:	Store at 4°C, do not freeze. Freezing R-Phycoerythrin conjugates will result in substantial loss of activity.
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: Neat; 10 ul labels 5x10 ⁵ cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences