

MHC Class I H-2Kk Mouse anti-Mouse Monoclonal (FITC) Antibody - LS-C72577 - LSBio

CatalogID:	LS-C72577
Target:	MHC Class I H-2Kk
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse
Antigen Type:	Cells
Immunogen:	A. AL mouse splenocytes.
Specificity:	Reacts with the H-2Kk MHC class I alloantigen. Cross-reactivity with splenocytes of SJL/Hsd mice has been observed by flow cytometric analysis. Does not react with other (e. g., b, d, q) haplotypes.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Complement-dependent cytotoxicity. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences